



# Gold Anchor™

## Product overview

Indicated use: The Gold Anchor™ marker is indicated for use to radiographically mark soft tissue for future therapeutic procedures.

The Gold Anchor™ needle is delivered pre-loaded with the Gold Anchor™ marker inside. It comes in a blister single pack, sterilized, ready for use.

Products are packed and sold in boxes of 10 units per product code.



## Minimally invasive insertion needles

The Gold Anchor™ marker comes preloaded in fine needles of three different diameters, 25G, 22G and 20G, with centimeter marking, adjustable depth limiter, echogenic surface and triple-sharpened gamma needle tip. These industry leading needle dimensions reduce the risk of complications during implantation of the marker.




## Markers

The Gold Anchor™ marker is available in two different diameters, 0.28 and 0.40 mm, and two different lengths, 10 and 20 mm, resulting in four different versions.

Depending on marker implantation technique, the markers can be implanted with a line shape or ball shape.

Choice of marker and marker shape depends on the user preferences regarding visibility, artifacts and dose perturbation.



## Product codes

Needle Ø	● 25G (0.5 mm)	● 22G (0.7 mm)				● 20G (0.9 mm)	
Length (cm)	15	8	15	20	25	20	35
<b>Marker</b>							
10	GA150-10			GA200-10	GA250-10		
20	GA150-20			GA200-20	GA250-20		
10-B		GA80-10-B	GA150X-10-B	GA200-10-B	GA250-10-B	GA200X-10-B	GA350-10-B
20-B		GA80-20-B	GA150X-20-B	GA200-20-B	GA250-20-B	GA200X-20-B	

Naslund Medical AB  
Åvägen 40B  
14130 Huddinge  
Sweden

Naslund Medical Inc.  
150 N Michigan Ave Ste 1950  
Chicago, IL 60601  
USA

www.GoldAnchorMarker.com  
Email: info@GoldAnchorMarker.com  
Order: order@GoldAnchorMarker.com